

is a 5cm (or larger or smaller) diameter removable disc to empty particle contents on the entry side and evacuate fluid from the entire reservoir, and to clean it. The opposite hole is a large diameter "female" entry port which accepts the extension tubing via an adaptor which importantly keeps the entry to the reservoir a large diameter and is centered from the edge of that side.

5. We claim the Khan Orotracheal Suction system can be used with the reservoir or hooked directly to standard commercially available suction containers.

Preliminary Amendment to Specification

Please place at the end of the specification. The Khan Orotracheal Suction system can be used with the reservoir or hooked directly to standard commercially available suction containers.